
# LightNeedle

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

\* Patients' sensation during and after laser needle versus metal needle treatment". van Amerongen KS, et al  
RESULTS: The common metal needle technique was well known by the patients in comparison to the laser needle method ( $p < 0.0001^{***}$ ). Laser needle acupuncture is a method which is painless ( $p < 0.0001^{***}$ ), energy inducing and relaxing ( $p = 0.0257^*$ ) which leads to a warming sensation ( $p = 0.0009^{***}$ ) during treatment.


\*\* Spectroscopy 16 (2002) 335-342, Seite 6  
Near-infrared spectroscopy for objectifying cerebral effects of needle and laser needle acupuncture  
Gerhard Litscher, Department of Biomedical Engineering and Research in Anesthesia and Critical Care, University of Graz, Auenbruggerplatz 29, A-8036 Graz, Austria  
Detlef Schikora, Department of Physics and Optoelectronics, University of Paderborn, D-33095 Paderborn, Germany

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Effect of 655-nm Low-Level Laser Therapy on

Exercise-Induced Skeletal Muscle Fatigue in Humans

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



Siehe auch;

A systematic review of low level laser therapy with location-specific doses for pain from joint disorders.

A systematic review of low level laser therapy with location-specific doses for pain from chronic joint disorders

Jan M Bjordal<sup>1</sup>, Christian Couppé<sup>2</sup>, Roberta T Chow<sup>3</sup>, Jan Tunér<sup>4</sup> and Elisabeth Anne Ljunggren<sup>1</sup> <sup>1</sup>University of Bergen, Norway <sup>2</sup>Lund University, Sweden <sup>3</sup>Private Medical Practice, Sydney <sup>4</sup>Private Dental Practice, Stockholm, Sweden

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## **Laser-Needle Therapy for Spontaneous Osteonecrosis of the Knee**

Winfried Banzer, M.D., Ph.D.,<sup>1</sup> Markus Hübscher, Ph.D.,<sup>1</sup> and Detlef Schikora, Ph.D.<sup>2</sup>

**Objective:** This case report describes the treatment of a 63-year-old patient with spontaneous osteonecrosis of the knee (SONK).

**Background Data:** SONK usually appears in the elderly patient without the typical risk factors for osteonecrosis. It is characterized by acute and sudden pain, mostly occurring at the medial side of the knee joint. Symptoms usually worsen with physical activity and improve with rest. Besides physical therapy, limited weight-bearing and the use of analgesics and nonsteroidal anti-inflammatory drugs, we propose low level laser therapy (LLLT) as a conservative treatment option.

**Methods:** LLLT was carried out using laser needles emitting radiation with wavelengths of 685 and 885 nm, and a power density of 17.8 W/cm<sup>2</sup>. Therapy sessions lasted 60 min and were performed daily over a period of 3 mo. The total irradiation dose emitted by 8 laser needles in 60 min of treatment was 1008 J. **Results:** Magnetic resonance imaging revealed distinct restitution

of the spongiosa edema 5 wk after treatment onset, and the final check-up at 35 wk demonstrated complete restoration of integrity. **Conclusion:** The present case report provides the first indication that laser-needle therapy may be a promising tool for complementary and alternative therapeutic intervention for those with SONK.

# Abstracts - Laserakupunktur

## **Patients' sensation during and after laserneedle versus metal needle treatment.**

van Amerongen KS, Kuhn A, Mueller M.

Department of Obstetrics and Gynaecology, Inselspital, Bern University Hospital, and University of Bern, Switzerland. k.staehler@spin.ch

**OBJECTIVES:** Aim of the study was to evaluate the patients' sensations during and after laserneedle versus metal needle acupuncture.

**STUDY DESIGN:** The prospective study was performed at the gynaecological outpatient department of a University Teaching Hospital of Bern, Switzerland. Thirty female patients per group were included in the study and randomized into laserneedle or metal needle group. All women visited the acupuncture out patient department because of gynaecological disorders. Age of the patients in the metal needle group was 38 years in median (range 18-73 years); mean age was 41+/-13.3. Age in the laserneedle group was 36 years in median (range 16-60 years) and mean age was 39.1+/-12.2.

Interventions were laserneedle acupuncture and metal needle acupuncture. Patients answered a questionnaire before, after the first treatment and prior to the second treatment. The questionnaires asked about the patients' knowledge of the various acupuncture methods and their health condition before treatment, their perception of pain, warmth, tiredness and relaxation during or after application of the needles or during or after the treatment. Statistics were performed by Graph Pad InStat 3 for windows.

**RESULTS:** The common metal needle technique was well known by the patients in comparison to the laserneedle method ( $p<0.0001^{***}$ ). Laserneedle acupuncture is a method which is painless ( $p<0.0001^{***}$ ), energy inducing and relaxing ( $p=0.0257^*$ ) which leads to a warming sensation ( $p=0.0009^{***}$ ) during treatment.

**CONCLUSION:** Both methods laserneedle and metal needle acupuncture are valuable methods in achieving relaxation and improvement of gynaecological symptoms. Laserneedle acupuncture is painless and easy to apply which is a valuable reason to support this technique in the future.

PMID: 18986749 [PubMed - indexed for MEDLINE]

# Laserblutbestrahlung

## **The application of intravenous laser irradiation of the blood for the combined treatment of psoriasis.**

Vopr Kurortol Fizioter Lech Fiz Kult. 2012 Jan-Feb;(1):24–8. Shakhova AS, Kulikov AG, Korsunskaja IM.

The objective of the present comparative study was to estimate the efficacy of the treatment of 104 patients presenting with psoriasis with the use of PUVA therapy and intravenous laser irradiation of the blood or the combination of the two methods. It was shown that the latter treatment caused a faster reduction of the inflammatory processes than the former whereas their joint application was especially efficacious in inducing regression of concomitant infiltrative events. Erythrocytes from the patients undergoing intravenous laser irradiation of the blood were characterized by the enhanced activity of antioxidative enzymes that are known to play an important role in the cellular protection from oxidants. At the same time, PUVA therapy caused no appreciable changes in the antioxidative status. Both treatments increased the levels of antiinflammatory cytokines and decreased those of proinflammatory ones which suggests their immunomodulatory effect that was especially pronounced in the case of combined therapy. It is concluded that the introduction of intravenous laser irradiation of the blood in the combined treatment of psoriatic patients exerts positive effect on compromised immunoregulation, stimulates the antioxidative system, and improves microcirculation.

PMID: 22693740[PubMed – indexed for MEDLINE]

## **Intravascular laser therapy (IVL) in pre-hypertension and hypertension conditions.**

S. Mokmeli MD, Sh. Bishe, Kh. Kahe, M. Shakhsh Laser Department of Milad Hospital. Hemmat Exp way , Chamran Exp. Way , Tehran , 119/ 1449614531, Iran, mokmeli@gmail.com This study was presented in In 7th international WALT congress, south Africa, 2008

Introduction: „Hypertension,“ is a condition that afflicts almost 1 billion people of worldwide and it is called the silent killer because it usually has no symptoms. Some people may find it out when they have trouble with their heart, brain, or kidneys. Intravenous laser therapy (IVL) is a technique of low level laser therapy that is pioneered in Russia and refers to blood irradiation through an intravenous needle or transcutaneous venous irradiation. IVL increases oxygen content, reduces the partial pressure of carbon dioxide and thrombocyte aggregation, stimulates fibrinolysis, and normalizes tissue metabolism... The aim of our study was evaluation the immediate effects of IVL in different grades of blood pressure.

Material and Method: Nowadays, since attention has been focused on systolic BP as a predictor of coronary and cerebrovascular disease, we allocated 125 patients according to the systolic blood pressure in 3 groups: 1-Normotensive (<120 mmHg n=50), 2-Pre-hypertensive (120–139 mmHg n=50), 3- hypertensive (stage I =140–159 mmHg, n=25). All the groups were conducted for 30 minutes IVL with a 630 nm, continues laser in 2.5 mW power at the end of intravenous fibre. Pulse rate, systolic, diastolic, and pulse pressures were measured before, after, and 15 minute after the IVL. All the results were analyzed by spss program.

Results: There was no statistically significant difference for pulse rate, systolic and diastolic blood pressure in normotensive group, however the significant difference was observed for pulse rate, systolic and diastolic blood pressure in Pre-hypertensive group as for systolic and diastolic blood pressure in hypertensive group(  $p < 0.005$ ).

Conclusion: IVL is an effective method for modifying factors to result a reduction in arterial pressure. It can be combined with anti hypertensive drugs in Pre-hypertensive and hypertensive patients as a modality of treatment; also it is a safe method in normotensive patients.

## **The development of method of intravenous laser irradiation of blood with green laser in patients with hyperlipidemia**

lin Lab Diagn. 2011 Aug;(8):15–7.

Poluéktova MV, Kharchenko IL, Kaplan MA, Sokol NI, Ershova LM, Borgul' OV, Chirkova TV, Vorob'eva OA.

The impact of intravenous laser irradiation of blood with green laser in patients with hyperlipidemia was investigated. The blood of patients was chosen as sample for analysis. The patients were divided in two groups: patients with atherosclerosis of various localization and patients with atherosclerosis associated with diabetes mellitus. The effectiveness of laser impact was evaluated according the blood biochemical indicators. The levels of crude cholesterol, triglycerides, low and very low density lipoproteins, apoproteins A and B, highly sensitive C-reactive

# Laserblutbestrahlung

protein, atherogenicity indicator, glucose content, uric acid content were determined before and after 1, 3 and 6 months after impact. The study results indicate the occurrence of hypolipidemic and hypoglycemic effects. PMID:22164411[PubMed – indexed for MEDLINE]

## **Effect of intravenous laser irradiation of the blood on plasma content of ceruloplasmin in patients with chronic pancreatitis.**

Bull Exp Biol Med. 2010 Nov;149(6):697–8. Burduli NM, Gutnova SK.

Department of Therapy, Advanced Training and Professional Re-Training Faculty, North Ossetian State Medical Academy, Federal Agency for Health Care and Social Development, Vladikavkaz, Russia.

### **Abstract**

Ceruloplasmin content in the blood plasma significantly increased in patients with chronic pancreatitis during exacerbation of the disease. Addition of intravenous laser irradiation of the blood to complex therapy of patients with chronic pancreatitis normalized ceruloplasmin content in the plasma.

PMID:21165421[PubMed – indexed for MEDLINE]

## **Intravenous laser irradiation of blood in the complex treatment for acute and progressive tuberculosis in adolescents**

Rusakova, L. I.; Dobkin, V. G.; Ovsyankina, E. S.; Firsova, V. A.; Bondarev, G. B.; Sidorov, V. A. Proc. SPIE Vol. 4422, p. 102–105, Low-Level Laser Therapy, Tatiana I. Solovieva; Ed.

Twenty-three out of forty-four adolescents affected by acute and progressive tuberculosis were administered a complex treatment based on the intravenous laser irradiation of blood. This treatment started 2 to 4 weeks after the beginning of chemotherapy. The efficiency of the treatment was assessed by clinical and laboratory criteria obtained after 1.5 to 2 months. It was found that the application of (lambda) equals 0.63 micrometers laser irradiation increased the therapeutic effect and stimulated the healing of the disease. Moreover, it enabled the achievement of a quiet course of tuberculosis with less pronounced residual changes in the lungs.

It has been studied the intravenous laser irradiation of blood (ILIB) influence with helium-neon laser (HNL) of 630 nm wavelength on some of lipid peroxidation (LPO) and antioxidant system (AOS) findings, aside-base status (ABS) and blood oxygen transport (BOT), state of dermal microhaemodynamics (MGD) in the intact rabbits and after modeling of local ischemia of brain (LIB). Depending on conditions of organism functioning (norm or brain ischaemia) ILIB has resulted in stimulating or normalizing effects on the whole metabolic and microhaemocirculation processes which had been studied during our investigation. It is discussed the mechanisms of pathogenetic directivity of ILIB influence in cerebral ischaemia

## **Intravenous laser irradiation of blood for the combined treatment of patients presenting with chronic sensorineural hearing loss**

Vestn Otorinolaringol. 2011;(2):43–5. [No authors listed] [Article in Russian]

A method for the treatment of chronic sensorineural hearing loss (CSNHL) is proposed that includes administration of trental followed by intravenous laser irradiation of blood (IVLBI). The study included 81 patients at the age from 20 to more than 60 years presenting with CSNHL. They were allocated to three groups; the patients in group 1 (n=32) were given trental intravenously followed by intravenous laser irradiation of blood, those of group 2 (n=24) were treated with IVLBI alone while patients of group 3 (n=25) received „traditional" treatment. Audiometric examination and rheoencephalography were carried out before and after therapy.

The hearing improved to 18–20 dB in group 1, to 10–15 dB in group 2, and to 10 dB in group 3. The improvement of rheoencephalographic characteristics was documented in the patients of all groups but was especially pronounced in group 1. It is concluded that the proposed method significantly increases the efficiency of treatment of chronic sensorineural hearing loss. PMID:21512485[PubMed – indexed for MEDLINE]

# Laserblubestrahlung

## **Antioxidant action and therapeutic efficacy of laser irradiation of blood in patients with ischemic heart disease**

Vopr Kurortol Fizioter Lech Fiz Kult. 2003 May-Jun;(3):22-5

Volotovskaia AV, Ulashchik VS, Filipovich VN.

Laser irradiation in therapeutic doses (gamma = 632.8 nm, 14 mW) has an antioxidant effect in blood irradiation in vitro as shown by activation of superoxide dismutase (SOD) which is a key enzyme of the antioxidant system (AOS) and suppression of lipid peroxidation. Adjuvant supravascular He-Ne laser irradiation of blood in combined therapy of 82 patients with ischemic heart disease (IHD) produces a positive trend in the clinical picture, hemostasis, lipid metabolism, blood SOD activity. Thereby, this method of laser hemotherapy is recommended for use in IHD patients. The dependence of the treatment results on the initial blood AOS necessitates consideration of AOS state in deciding on laser therapy in this group of patients.

PMID: 12852010 [PubMed - indexed for MEDLINE]

## **A possible mechanism of the hypotensive effect of laser irradiation in patients with ischemic heart disease with arterial hypertension**

Vopr Kurortol Fizioter Lech Fiz Kult. 2001 Sep-Oct;(5):17-20.

Vasil'ev AP, Strel'tsova NN.

93 patients with ischemic heart disease, 44(47.3%) of them with arterial hypertension, were exposed to laser irradiation (LI). LI was accompanied with lowering of arterial pressure, more prominent in hypertensive patients. Mechanism of the hypotensive effect of laser therapy operates largely through stabilization of the lipid bilayer of the cell membrane demonstrated on the model of erythrocyte. PMID: 11785330 [PubMed - indexed for MEDLINE]

## **Efficiency of low-intensity laser radiation in essential hypertension**

Velizhanina IA, Gapon LI, Shabalina MS, Kamalova NN. Klin Med (Mosk). 2001;79(1):41-4

In a placebo-controlled study an antihypertensive activity of low-intensive laser radiation (LILR) was evaluated in 52 males with essential hypertension stage I. The placebo group consisted of 14 matched patients. LILR was used as monotherapy of 10 daily procedures. This treatment significantly lowered systolic, diastolic and mean arterial pressure. Moreover, diastolic arterial pressure did not rise high at submaximal bicycle exercise. Total peripheral vascular resistance also decreased. A good hypotensive effect was achieved in 90.4% cases. Thus, LILR is a highly effective treatment in essential hypertension stage I.

Publication Types: Clinical Trial PMID: 11234266 [PubMed - indexed for MEDLINE]

## **The blood antioxidant system in patients with ischemic heart disease undergoing laser therapy**

Vopr Kurortol Fizioter Lech Fiz Kult. 1993 Nov-Dec;(6):4-7.

Zubkova SM, Sorokina EI, Kenevich NA, Tupitsyna Iulu, Minenkov AA.

Catalyse activity and lipid peroxidation were studied in the sera of patients with coronary heart disease who had been exposed to red and infrared low-energy laser radiation (LLR) applied to the skin projections of the heart and its reflexogenic areas. It was found that the optimal normalizing effect occurred with the red band of LLR applied in the rehabilitative period to patients with prior myocardial infarction concurrent with hypertensive disease and stable angina of functional classes I and II. In functional class III stable angina, infrared LLR had advantages over the red on when they acted on the blood antioxidative system. PMID: 8154120

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REIMERS & JANSSEN GmbH